

ALARA in Handling of Beta-Emitters - Measurement Techniques and Optimisation

Ilona Barth, Arndt Rimpler

Bundesamt für Strahlenschutz, SG 2.2, Köpenicker Allee 120-130, 10318 Berlin



13th EAN Workshop, 7-10 June 2011, Oscarborg Fortress (Norway)



Introduction

- Increasing use of beta emitters in nuclear medicine diagnostic and therapy
- Necessity to handle high activities at small distances to the skin
→ increased probability of high skin exposure (>500 mSv/a)!
- Lacking individual awareness of high skin doses to the staff resulting from absent, inadequate or inaccurately placed extremity dosimeters in routine monitoring
→ low radiation protection standard

Introduction

Therapies under review:

- Radiosynoviorthesis (RSO) using Er-169, Re-186, Y-90
- Radioimmunotherapy (RIT) using Y-90-Zevalin[®]
- Radiopeptidtherapy (PRRT) using Y-90-DOTATOC
- Intravascular Brachytherapy (IVB) using Re-188
- Palliative Therapy using Sm-153

Methods and Materials



Thin-layer thermoluminescence detectors (TLD)

Material:

LiF:Mg,Cu,P type: MCP-Ns™
Area mass of the sensitive layer
<10 mg/cm²

Calibration:

Beta-Secondary Standard facility
BSS 2 with a Sr-90/Y-90-source

Evaluation:

Harshaw 3500 type reader

Results / Radiosynoviorthesis with Y-90

Measurements performed in:

- 11** Doctor's surgeries
- 13** Radiologists / technicians
- 18** Doctors
- 210** Patients
- 45** GBq Y-90 administered




Results / Radiosynoviorthesis with Y-90

Influence of the shielding on skin exposure during syringe dispensing

without shielding

wrong !

$\approx 50 \mu\text{Sv}/\text{MBq}$



with shielding for vial, syringe
and upper needle


$\approx 2 \mu\text{Sv}/\text{MBq}$




Results / Radiosynoviorthesis with Y-90

Influence of the Makrolonring on the skin exposure

without Makrolonring




with Makrolonring



RSO with Y-90: local normalized skin dose

Results / RSO, Administration

Causing very high exposure



Direct contact of the fingers to the needle: **~ 2 mSv/s!**

Results / European ORAMED-Project

SCOPE


Radiotherapies with Y-90

Scope of measurements in preparation (P) and administration (A) of Y-90 in NM therapies

Therapy	Procedures	Countries (B,Ch,D,E, F,I)	Number of		Data sets
			Hospitals	Workers	
PRRT <i>Dotatoc</i>	P	3	3	5	16
	A	3	3	7	17
RIT <i>Zevalin</i> ®	P	6	1	20	49
	A	6	6	22	45
			1		
			5		




RESULTS / Classification of workers for RIT with Y-90/Zevalin®, Preparation




Results / Zevalin® Preparation

RESULTS

Causing very high exposure during Y-90/Zevalin® Preparation





RESULTS/ Classification of workers for RIT, Administration



Results / Zevalin®

Causation of very high exposure



Results / Zevalin[®], Administration

Any direct contact to vessels with activity must and can be avoided!!!

Reduction of exposure more than one order of magnitude:





Results / ORAMED

Averaged maximum skin dose per preparation or administration in nuclear medicine therapies

Therapy	Procedure P Preparation A Administration	Worker	Max. Skin dose [mSv]			
			Mean	Median	Min	Max
⁹⁰Y - RIT	P	15	16.5	14.2	1.8	65.9
	A	19	4.8	2.9	1.0	11.9
⁹⁰Y-PRRT	P	5	21.6	11.3	1.0	76.2
	A	7	10.4	8.2	2.2	26.9

Mean activity per procedure: RIT-P: 1.5 GBq RIT-A: 1.0 GBq
PRRT-P: 10.3 GBq PRRT-A: 5.5 GBq

Results / RIT



Results / ORAMED

RESULTS

Why are/were the high skin doses at the fingers not detected in routine monitoring???

Possible reasons:



- Inadequate dosimeter?
- Inadequate position where the dosimeters are worn?

Answer:

- Ringdosimeter with thin layer TLD of [LiF: Mg,Cu.P] are appropriate to measure correct beta and photon radiation dose
- The best position of routine ringdosimeter should be established in ORAMED studies.

Results / ORAMED

Ratio of maximum dose on both hands to dose on dosemeter position



13th EAN Workshop, 7-10 June 2011, Oscarborg Fortress (Norway)



Results / ORAMED

- The best position for the routine ringdosemeter:
the base of the index finger of the **non-dominant** hand!
- The ratio of the maximum dose of both hands to the dose of the dosemeter on the base of the index finger of the non-dominant hand is about **6 for all procedures in nuclear medicine** (therapies using Y-90 and diagnostics using Tc-99^m and F-18). Generally, a worker performs more than one procedure!



Conclusions

- Local skin dose can exceed the limit (500 mSv/a) by numerous workers in hospitals where radiation protection standard is low
- There is a high potential to decrease exposures by simple means ► **shielding and tools**
- The **awareness** of the necessity to improve the radiation protection standard must be increased
- Adequate skin dose monitoring is urgently needed in nuclear medicine
- Ring dosimeters should be worn **on the base of the index finger of the non-dominant hand**
- Even if the dosimeter is worn on the base of the index finger, the skin dose maximum is considerably **underestimated by a factor of about 6**
- The ringfinger base of the dominant hand, where ring dosimeters are often worn, is not a proper position in routine monitoring and wrist dosimeters are not appropriate for monitoring nuclear medicine staff



Acknowledgement

Parts of the research leading to these results has received funding from the European Atomic Energy Community's Seventh Framework Programme (FP7/2007-2011) under grant agreement n° 211361.

Thanks to all members of medical staff and hospitals supporting the research.

Thank you for your attention!

